

## GUI Manual

*Load .pos file*

File > Load... > .pos

Cmd + L

*View undo/redo history*

View > Tools

Cmd + T

Select History tab

*Ion suggest*

View > Tools

Cmd + T

Select Peak ID tab

*Add ion(s)*

In the Peak ID tab, select ions (Ctrl + click) > Add

OR

In the spectrum, double-click one suggest line > press 'a'

*Delete ion*

In the table, select ions (Ctrl + click) > right-click > Delete

Del

*Add/edit range*

In the table, double-click the start and end values to edit

OR

In the spectrum, double-click one suggest line > press 's' to enable the span

> draw span > press 'c' to commit

*Export to .rng file*

Edit > Export... > .rng

Cmd + R

*Export to. json file*

Edit > Export... > .json

Cmd + J